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
- Cesar Espiritu, MD, vice-chairman of the Manila Doctors Hospital ophthalmology department

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
Lab testing conditions
 AMO* set at 95 cm bottle height
 CENTURION® set at 50 mmHg

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Lab testing conditions
 B&L* setting 70 cm bottle height with 70 mmHg added
 CENTURION® set at 50 mmHg

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


Group difference (95% confidence interval):
-2.79 (-3.44 to -2.13) percent-seconds;
data reflect least squares mean \pm standard error.

Group difference (95% confidence interval):
-6.12 (-9.82 to -2.43) mL; data reflect least squares mean \pm standard error.

Group difference (95% confidence interval):
-15.78 (-26.49 to -5.07) seconds; data reflect least squares mean \pm standard error.

Centurion[®] Configuration: Centurion[®] Vision System, 45° Balanced Tip with Intrepid[®] Ultra Sleeve

INFINITI[®] Configuration: INFINITI[®] Vision System, 45° Mini Flared Kelman with Ultra Sleeve

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CENTURION[®] Vision System with ACTIVE SENTRY[®] Handpiece:

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*As compared to the INFINITI[®] Vision System, bottle gravity system.

1. Sharif-Kashani P, Fanney D, Injev V. Comparison of occlusion break responses and vacuum rise times of phacoemulsification systems. *BMC Ophthalmol.* 2014;14:96.
2. Alcon data on file. 2017
3. Nicolini CM, Dimalanta R, Miller K. Experimental anterior chamber maintenance in active versus passive phacoemulsification fluidics systems. *J Cataract Refract Surg.* 2016;42(1):157-162.
4. Alcon data on file. 2017
5. Khokhar S, Aron N, Sen S, Pillay G, Agarwal E. Effect of balanced phacoemulsification tip on the outcomes of torsional phacoemulsification using an active-fluidics system. *J Cataract Refract Surg.* 2017;43(1):22-28.
6. Solomon K, Lorente R, Cionni R, Fanney D. Prospective, randomized clinical study using a new phaco system with intraocular system target pressure control. ASCRS-ASOA Symposium and Congress; April 25-29, 2014; Boston, MA.
7. Zacharias J. Comparative motion profile characterization of the mini flared and balanced phacoemulsification tips. ESCRS; September 5-9, 2015; Barcelona, Spain.
8. Alcon data on file. 2017
9. Alcon data on file. 2017
10. Aravena C, Dyk D, Thorne A, Fanney D, Miller K. Percent aqueous volume loss associated with post occlusion break surge in 4 phacoemulsification systems. ASCRS-ASOA Symposium and Congress; May 5-10, 2016; New Orleans, LA.
11. Boukhny M, Sorensen G, Gordon R. A novel phacoemulsification system utilizing feedback based IOP target control. ASCRS-ASOA Symposium and Congress; April 25-29, 2014; Boston, MA.
12. Zacharias J. Comparative thermal characterization of phacoemulsification probes operated in elliptical, torsional and longitudinal ultrasound modalities. ASCRS-ASOA Symposium and Congress; April 25-29, 2014; Boston, MA.
13. Vasavada AR, et al. Comparison of torsional and microburst longitudinal phacoemulsification: a prospective, randomized, masked clinical trial. *Ophthalmic Surg Lasers Imaging.* 2010;41(1):109-114.
14. Allen D, Habib M, Steel D. Final incision size after implantation of a hydrophobic acrylic aspheric intraocular lens: new motorized injector versus standard manual injector. *J Cataract Refract Surg.* 2012;38(2):249-255.
15. Johansson C. Comparison of motorized IOL insertion to traditional manual IOL delivery. ASCRS-ASOA Symposium and Congress; March 25-29, 2011; San Diego, CA.

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